



## Anterior Restoration Checklist

- Detailed lab slip including:**
  - Teeth to be restored with type of restoration
  - Shade selected
  - Stump shade for all-ceramic restorations
  - Material choice
    - for PFMs select the class of metal, such as noble, high noble, etc.
  - Patient goals
- Impression of maxillary and mandibular arches**
- Pre-op models**
- Bite record**
- Model or impression of approved provisionals**, when necessary
- Photographs**, suggested views:
  1. Full face with natural smile
  2. 1:2 full smile
  3. 1:2 retracted, teeth slightly parted, frontal, left and right views
  4. 1:1 close-up showing any surface texture, incisal details, teeth parted, slightly underexposed.
  5. Shade tab held up to the incisal edge and in the same plane as the tooth, with the shade tab label clearly visible.
  6. Stump shade tab held next to prepped tooth.
  7. Any other views or videos which may be important to the successful completion of the case.

Photographs can be sent in several ways:

1. Compressed and attached to an email.
2. Burned onto a CD or loaded onto a memory stick. These will be returned with the completed case.
3. Post photos to a share photo site, such as Dropbox.